

# YDM-PRO AMALGAMATOR

YDM-pro is a new type mixer for all mixed capsules including amalgam capsules. The universal high-speed amalgamator is controlled by a microprocessor, which ensures a precise, controllable and consistent mixing. With a safety cover, YDM-pro will deactivate automatically when you open the front cover to access the capsule.

## FEATURES:

- Silicone rubber buttons for comfortable touch
- LED display
- Two speed selections
- Noiseless and stable operation
- Easy operation and sturdy capsule holder
- Sealed mixing area with a safety cover
- Ergonomic design



## TECHNICAL DATA:

- Mixing frequency: 4500 and 4000 oscillations per sec +/- 5%
- Mixing time: Adjustable from 1 to 99 seconds
- Operating temperature: 5 to 40°C
- Voltage / Current: 200 TO 240 Volts, 50 Hz / 30 Watts
- Dimensions: 205L\*195W\*115\*(without power cord)
- Weight: 2.2Kg
- Amplitude: 15mm
- Functioning: intermittent