

E103 Dental Autoclave



Automatic

Chamber Size	Ø250XD350mm
Chamber Capacity	18L
Autoclave Outer Size (L*W*H)CM	56*39*36
Package Size(L*W*H)CM	65*46*41 CM / W:24KGS G.W: 26.5KGS
Sterilize Temp.	122°C & 132°C
Sterilize Pressure (Mpa)	0.12Mpa & 0.23Mpa
Sterilize Time	20Min & 6Min
Power Consumption	900W
Material	Stainless Steel SUS#304.
Current -110V	8 A
Current -220V	4 A
Auto Water Supply Auto Drainage	Yes
Safety configuration	*Safety Valve: if inner pressure over standard, vent automatically. *Overheat Protect Device: when inner temp overheat, power off automatically
Water Store Tank	Outer
Water Feeding Pump	Yes
Temp. Disply Window	LED display window(no operation symple) and pressure gage
Sterilize Principle	Uses the vapor power emissions cold air / 2 time