

PORT-X III

Portable Dental X-ray System

- Excellent Image Quality
- Elegant Design
- Perfect Output for Film, Digital Sensor : 60kV, 2mA
- Easy Operation



Compact Outside, Smart Inside

New paradigm of Intra-Oral X-ray system is now at hand.
PORT-X III makes your dental clinic full of ease and comfort.

Digital X-ray system maker, **GENORAY Co.,Ltd.**

#512 Byucksan Technopia 434-6, Sangdaewon 1-Dong, Jungwon-Gu, Seongnam-City, Gyeonggi-Do, KOREA
TEL : +82-31-737-8020 FAX : +82-31-737-8025 E-mail : genoray@genoray.com <http://www.genoray.com>



PORT-X III

Portable Dental X-ray System

Features

- Light and small body with rechargeable which enables you to shoot easily and to keep everywhere.
- Small does of radiation keeps patients and dentists away from radiation.
- High performance Tube and HFG make sharp and clear images.
- Use with digital sensor as well as chemical film.



Components of PORT-X III


Main Unit
 - X-ray tube assembly / Beam limiting Device
 / Battery

Accessories
 - Adaptor / Staying Necklace

Technical Specification & Dimension

X-ray Generator	Type	High Frequency Inverter
	Power Output	180W
X-ray Installation	Tube Voltage	60 kV
	Tube Current	2mA
	Heat Capacity	8.5 KHU
X-ray Tube	Anode Type	Stationary
	Anode Angle	20°
	Focal Spot	0.8mm
X-ray Control	Exposure Time Set	46 Steps (0.01 ~ 2.0 sec)
Power Requirement		22.2 VDC
Weight	Main Body	2.35 kg
	Accessory	0.6 kg
Dimension	Main Body	197 x 147 x 145 (mm)
	Charger	100 x 65 x 35 (mm)

• Please do not hesitate to contact us, if you need the detailed specification of our products.



Graphic LCD
 Users can view the operation through graphic LCD with very high quality.

Operation Panel
 Users can make the most optimized x-ray control suitable to patients and the best images just through simple key operation.

Digital X-ray system maker, **GENORAY Co.,Ltd.**

#512 Byucksan Technopia 434-6, Sangdaewon 1-Dong, Jungwon-Gu, Seongnam-City, Gyeonggi-Do, KOREA
 TEL : +82-31-737-8020 FAX : +82-31-737-8025 E-mail : genoray@genoray.com <http://www.genoray.com>

