

Port-XIV®

Portable Dental X-ray System

1.5 kg
Light weight

0.4 mm
Focal spot

70 kv
Tube voltage



Superior image quality
Portable camera type
Easy to control UI
Compact size

Compact and sharp portable X-ray

- Ergonomic design for one handed operation
- Lightweight and compact due to optimal tank design
- Image viewer with clear IPS LCD Screen **VIPS**
- Winner of Red Dot Award for Design 2017


Easy to use portable X-ray

- Straightforward process for image capture with intuitive icons and 5-ways button
- High-capacity battery capable of lasting all day without the need to recharge Li-polymer 22.2VDC, 1,000mAh

Superior image quality

- Crystal clear images from 0.4mm Focal spot and 70kv X-ray Tube

▶ Effect of Focal Spot Size



Smart Portable X-ray



◦ 5-ways button

Setting X-ray exposure time and selecting teeth or menu and the other options by moving up, down, left and right.

◦ LCD Screen

Set and Check the exposure options through LCD screen.

◦ Hand strap

A hand-strap makes stable exposure without trembling.

◦ Key

Quickly and easily select exposure mode and menu setting.



GENORAY

T +82-31-5178-5500 E genoray@genoray.com W www.genoray.com

CE
2460

* This X-ray unit may be dangerous to patient and operator unless safe exposure factors and operating instructions are observed.

© GENORAY CO.,LTD. Ver. 2023-12