

# The New Worldwide Standard for Equipment Support and Radiation Protection Systems – PORTEGRA2



**PORTEGRA2**  
designed by MAVIG

PATENTED by **MAVIG**

- Ceiling anchor
- Brake mechanism for track-guided systems
- The MAVIG spring arm
- Internal safety spring
- The safety notch

**SAFETY FEATURES**

- Single piece, formed-steel ceiling column
- Strengthened wall thickness at the mounting post
- Maximum load capacity of 18 kg (40 lbs) on each mounting post
- The steel extension/spring arm ensures high stability against collisions
- UL listing and CE certification
- User-friendly, smooth edged design
- Allows for different loads levels
- Easy spring-arm tension adjustment
- Life Cycle Testing