



EMBRACE EXCELLENCE



THERMOSEALER – E301

- 1. Ergonomic design
- 2. Desktop
- 3. Compact size
- 4. Practical and easy to use
- 5. Quick sealing
- 6. CE certificate
- 7. Hidden Blade more safe and reliable



Packet Tray		Optional
Overall dimension WXHxD	cm	51X45X19
Electrical Voltage and frequency	V-HZ	110-230/50(single-phase AC)
Absorb maximal power	W	100
Seal width	mm	12
Seal length	mm	320
Seal duration	S	3
Weight	Kg	7.5