

M203 Endo Motor with apex locator

Feature:

R-Smart Plus New Endo Motor
With Apex Locator Function
German Motor
Large Colorful OLED Screen
Four Working Modes
Six Functions

Technical:

Technical specifications:

Battery: 7.2V/1400mAh

Adapter: 100V~240V 0.8A 50Hz/60Hz

Revolving speed: 200~600rpm

Torsion: 0.6~4.0N.cm

Buzzer alert: The buzzer will alert when the file is less than 2mm closed to the apex

Dimensions: 250mm*210mm*140mm

Package weight: 1.45kg

