

N103 LED Ultrasonic Scaler with Water Bottle



Features:

- *Automatic Frequency Tracking
- *Two Water Supply: Auto-water supply OR external
- *Control: Wireless foot switch
- *Autoclavable detachable handpiece
- *Imported long life liquid pump
- *Military use power tube inside, more reliable
- *Titanium alloy acoustic horn, anti-fatigue
- *Soft anti-wear cable
- *Handpiece: Illuminate detachable handpiece
- *Function: Scaling, Perio, Endo irrigation
- *Handpiece Detachable handpiece
- *Functions Scaling&Endo&Perio
- *Control mode Wireless control
- *Specifications:
- *Foot switch battery 1.5v*2
- *Receiving sensitivity -1140b
- *Water pressure 0.01MPa-0.05MPa
- *Input AV100-240V 47-63Hz
- *Power 30VA-48VA
- *Main unit input DC30V 1.5A
- *Output half-excursion <math>< 2N</math>
- *Vibration frequency of the tip 25~31KHz Primary
- *Vibration excursion of the tip $\leq 100\mu\text{m}$
- *Ultrasonic output power 3W-20W
- *Weight of main unit:1.10KG
- *Gross Weight:2.5KG
- *Dimensions:330mm*155mm*155mm
- *Packing Size:340*240*180mm

PRODUCT ADVANTAGE

Wireless control



- ▶ no connection between scaler and foot switch
- ▶ easier for dentist to operate
- ▶ makes clinic space clearer and safer
- ▶ anti-interfere wireless technology
- ▶ need only 2 AA battery for one year using



Intelligent vibration technology



QUALITY MATERIALS

LED DETACHABLE HANDPIECE



pure white light, clear eyesight

autoclavable under high temperature and high pressure

small volume, light weight, easy operation

six bulbs inside LED lamp, two LED lamps for replacement

PRODUCT DISPLAY



