

GETIDY DYS – MOBILE X – RAY UNIT



DYS Series Dental X-Ray Unit

> Mobile Unit DYS-M

■ ■ ■ colour

Fashion Design
Stable and Strong Support Base
Strong column
Stretchable arm
Multi-function timer
Easy-operational handset

Specification:

Power Voltage: AC 220V \pm 10%
Frequency: 50/60Hz
Electric Plug: Three-core(international/
Europeanstandard)
Ultimate Capacity: 900VA
Electric Current: 4A
Fuse: 6.3A
Focus: 0.8mm
Rated High-Voltage: 70KV \pm 10%
Tube-Head Current: 7mA \pm 15%
Anode Angle: 19°
Load Cycle: 1/60
Half Value Layer: 1.6mmAL under 70kv
Inherent Filtration: \geq 2.1mmAl
Radiation Leakage: <0.07mGy/h
Exposure Time: 0.06-2.00s

> Package1: 695X385X475mm-25.3Kg
> Package2: 710X540X100mm-22.0Kg
> Package3: 1505X140X100mm-9.2Kg



> Fashion and Flexible
High-voltage generator



> Flexible Joint of Arm
Structure



> Fashion and Comfortable
Patient Seat



WE LOVE WHAT WE DO

