



GETIDY DYS – WALL MOUNTED X – RAY UNIT

DYS Series Dental X-Ray Unit



Installation Space

A=420mm (Standard) or 485mm (Optional)
 B=1055mm
 C=2230mm
 D=880mm



> Pack Size: 910X385X390mm-33.5Kg

Specification:

Power Voltage: AC 220V \pm 10%	Anode Angle: 19°
Frequency: 50/60Hz	Load Cycle: 1/60
Electric Plug: Three-core(international/ Europeanstandard)	Half Value Layer: 1.6mmALunder 70kv
Ultimate Capacity: 900VA	Inherent Filtration: \geq 2.1mmAl
Electric Current: 4A	Radiation Leakage: <0.07mGy/h
Fuse: 6.3A	Exposure Time: 0.06-2.00s
Focus: 0.8mm	
Rated High-Voltage: 70KV \pm 10%	
Tube-Head Current: 7mA \pm 15%	



WE LOVE WHAT WE DO

