



**INDEPENDENT X-RAY  
QUALITY ASSURANCE**

Use with

RTI Scatter Probe  Piranha  Cobia



**Ocean Next™ version 2.0** general specifications

## The world-leading software for X-ray QA & testing

Since its launch in 2010, Ocean™ has been the world-leading software for X-ray QA and testing, addressing customers' needs from Quick, for daily checks, to the advanced QA with workflow, automation, and traceability in Professional - in easy steps.

The Ocean Next™ software v.2.0 includes upgraded and new automated features, integration with myRTI, and much more! If you want a swift, easy-to-use application for routine controls, or a customized application with workflow, automatic tests, and traceability, Ocean continues to provide you with both options.

Ocean's recognition and efficiency make it central to our X-ray solutions, and it is compliant with our multimodality Piranha and Cobia meters and the RTI Scatter Probe combined leakage and scatter detector.

You can plan the measurements in advance, create checklists, add information as a pop-up window for a specific exposure, and include instructions to simplify the work for you and your co-workers using a streamlined user interface.

Which Ocean Next™ license suits you best? Choose from three different license levels Quick, Advantage, or Professional.

### Key features

- Customized application with workflow, automatic tests, and traceability.
- Updated, streamlined user interface for recognition and efficiency.
- Message panel for RTI news and more visible, automated calibration reminders.
- Calibration notices available in reports.
- Simplified templates and measurements structure.
- Runs on both computers and tablets - choose the display you want.
- RTI Academy training courses available.

### Which Ocean Next™ suits you best?

Built-in display function	Quick/Advantage/Professional
Save function	Quick/Advantage/Professional
Printing	Quick/Advantage/Professional
Basic Excel® connection	Quick/Advantage/Professional
Advanced Excel® connection	Advantage/Professional
User-defined tests & checklists	Advantage/Professional
Analysis with pass/fail limits	Advantage/Professional
Reporting	Advantage/Professional
CT Dose Profiler support	Advantage/Professional
Single-page templates	Advantage/Professional
Multi-page templates	Professional
Site database (including location, room, equipment, etc.)	Professional
Scheduling	Professional

**Still not sure? View the Ocean Next software page on the RTI website.**

RTI Group Headquarters  
Flöjelbergsgatan 8C  
SE-431 37 Mölndal, SWEDEN

Phone: +46 (0) 31 746 36 00  
sales@rtigroup.com  
www.rtigroup.com

RTI Group North America  
33 Jacksonville Road, Bldg. 1  
Towaco, NJ 07082, USA


Phone: +1 800-222-7537  
sales.us@rtigroup.com  
www.rtigroup.com



INDEPENDENT X-RAY  
QUALITY ASSURANCE

Use with

RTI Scatter Probe Piranha Cobia



### Ocean Next™ version 2.0 general specifications

## myRTI

myRTI is a free portal for all Ocean Next™ software version 2.0 users, providing many benefits for X-ray Quality Assurance, Quality Control testing, and measurements. Register with RTI today to access all the myRTI Key features and more!

### Key features to simple ownership

#### Logs & Calibration Certificates

Manage your RTI devices, access calibration certificates, and keep track of all your measurement data acquired using Ocean Next™. Accessible through myRTI from a PC anywhere!

#### Meter info




Manage your RTI devices using meter info. Track devices used with Ocean Next™, reminders on calibration due dates, and access to your calibration certificates.

#### Single sign-on account

Your myRTI account is a single sign-on account giving you access to Ocean Next™, the myRTI customer portal (meter info & logs), plus access to RTI Support and RTI Academy online training.

#### Support

One of our goals is to share our deep knowledge of best practices within the field of X-ray QA and testing – to ultimately protect patients and staff in an ever more complex operational environment.



RTI Group Headquarters  
Flöjelbergsgatan 8C  
SE-431 37 Mölndal, SWEDEN

Phone: +46 (0) 31 746 36 00  
sales@rtigroup.com  
www.rtigroup.com

RTI Group North America  
33 Jacksonville Road, Bldg. 1  
Towaco, NJ 07082, USA

Phone: +1 800-222-7537  
sales.us@rtigroup.com  
www.rtigroup.com